

Nicare System



Nicare special points

1. Multiple Digital System

- Combines all the advanced technology obesity treatment equipment
- The original system operated by WINDOWS
- Ability to store start-up parameters
 - Customer information and treatment mode
- Built in DSP, easy to change frequency, power, and range type
- Combination therapy(available with medium frequency and ultrasound)
- Both Dual channel

2. Innovative design

- Touch screen for easy selection of ultrasound and medium frequency controls
- Large multi-segment LCD display
- Easily adjustable flip-up LCD screen
- One touch movable keyboard operation simplifies ready to start therapy
- Support cord reel, for easy keeping, remote control power system
- Ultra slim design for ease of storage
- Support MP3function



Combination Obesity Treatment Equipment

Nicare Treatment



Medium Frequency



Ultrasound

Nicare Feature

MEDIUM FREQUENCY

- Pad type head with 2Khz, 4Khz, 8Khz frequencies
- Separate control for each pad to adjust power and therapy time
- Each frequency supports nine auto waveform mode types
- Supports two therapy types: "resolve fat" mode, "strengthen the muscles" mode
- Shipped with 12 sets of therapy pads


ULTRASOUND

- Sound head with 30 KHz
- Dual channel
- Output of power set up by **Continuance mode**/**User setup mode** type
- Includes 12 Sound heads, which can provide all over body therapy
- Support 12 channels

Nicare Accessories



Medium frequency acc



Ultrasound acc

